

Introduction

This Quick Reference Guide supplements the computer manuals included on the Getting Started CD included with the **Optimist MMX-3**.



Switch Inputs

There are two 3.5mm switch jacks on the back edges of the case. These inputs emulate the F1 (left rear) and F2 (right rear) keys. If your communication program requires other keyboard keys as switch inputs a different interface will be necessary, like a USB port or Bluetooth connection. Use normally open switches with these switch inputs.

Mounting Bracket

The optional clamp-on mounting bracket makes for easy attachment to wheelchair mounting systems, and holds the **Optimist MMX-3** safe and secure.

Resources

- >The Getting Started CD contains the **Optimist MMX-3** User's Manual.
- >For all **Optimist MMX-3C** devices with **The Grid**, the Getting Started CD also contains the Reference Manual, Getting Started, and the Training Cards for The Grid 2 program.
- >Web Site Support Downloads can be found at www.zygo-usa.com/Support/.
- >For additional support, contact ZYGO Customer Service

ZYGO-USA
48834 Kato Road, Suite 101-A
Fremont, CA 94538 U.S.A.
Web: www.zygo-usa.com

U.S.A./Canada: (800) 234-6006
Tel: +1 (510) 249-9660
Fax: +1 (510) 770-4930
E-mail: support@zygo-usa.com